THE FIRST AND ONLY

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- [1] Anderson RR, Parrish JA. Selective photothermolysis: precise microsurgery by selective absorption of pulsed radiation. *Science*. 1983;220(4596):524-527.
- [2] Freedman JR, Greene RM, Green JB. Histologic effects of resurfacing lasers. *Facial Plas Surg*. 2014;30(1):40-48.
- [3] Vogel A, Venugopalan V. Mechanisms of pulsed laser ablation of biological tissues. *Chem Rev*. 2003;103(2):577-644.
- [4] Tibor J, Kurtz R, Raksi F, Suarez C, Horvath C, Spooner G. The femtosecond blade: applications in corneal surgery. *Optics Photonics News*. 2002;13(1):24-29.

- [5] Girard B, Cloutier M, Wilson DJ, Clokie CM, Miller RJ, Wilson BC. Microtomographic analysis of healing of femtosecond laser bone calvarial wounds compared to mechanical instruments in mice with and without application of BMP-7. *Lasers Surg Med*. 2007;39(5):458-467.
- [6] Nagy ZZ. New technology update: femtosecond laser in cataract surgery. *Clin Ophthalmol*. 2014;8:1157-1167.
- [7] Dingus, RS, Scammon RJ. Grüneisen-stress-induced ablation of biological tissue. *Proc SPIE. Laser-Tissue Interaction II*. 1991;1427:45-54.
- [8] Paltauf G, Dyer PE. Photomechanical processes and effects in ablation. *Chem Rev*. 2003;103(2):487-518.
- [9] Ho DD, London R, Zimmerman GB, Young DA. Laser-tattoo removal—a study of the mechanism and the optimal treatment strategy via computer simulations. *Lasers Surg Med*. 2002;30(5):389-397.
- [10] Chen H, Diebold G. Chemical generation of acoustic waves: a giant photoacoustic effect. *Science*. 1995;270(5238):963-966.
- [11] Ross V, Naseef G, Lin G, et al. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. *Arch Dermatol*. 1998;134(2):167-171.
- [12] Herd RM, Alora MB, Smoller B, Arndt KA, Dover JS. A clinical and histologic prospective controlled comparative study of the picosecond titanium:sapphire (795 nm) laser versus the Q-switched alexandrite (752 nm) laser for removing tattoo pigment. *J Am Acad Dermatol*. 1999;40(4):603-606.
- [13] Kilmer SL, R. E. Fitzpatrick, L. B. Da Silva, R. Marshall, R. Ghiselli. "Picosecond and femtosecond laser treatment of tattoo ink." *Lasers Surg Med*. 1996;19(Suppl 8):36.
- [14] Izkison L, Farinelli W, Sakamoto F, Tannous Z, Anderson RR. Safety and effectiveness of black tattoo clearance in a pig model after a single treatment with a novel 758 nm 500 picosecond laser: a pilot study. *Lasers Surg Med*. 2010;42(7):640-646.
- [15] Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: a prospective trial. *Arch Dermatol*. 2012;148(12):1360-1363.
- [16] Brauer JA, Reddy KK, Anolik R, et al. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. *Arch Dermatol*. 2012;148(7):820-823.
- [17] Franjic K, Cowan ML, Kraemer D, Miller RJD. Laser selective cutting of biological tissues by impulsive heat deposition through ultrafast vibrational excitations. *Opt Express*. 2009;17:22937-22959.
- [18] Amini-Nik S, Kraemer D, Cowan ML, et al. Ultrafast mid-IR laser scalpel: protein signals of the fundamental limits to minimally invasive surgery. *PLoS One*. 2010;5(9).
- [19] Böttcher A, Clauditz TS, Knecht R, et al. A novel tool in laryngeal surgery: preliminary results of the picosecond infrared laser. *Laryngoscope*. 2013;123(11):2770-2775.
- [20] Jowett N, Wollmer W, Mlynarek AM, et al. Heat generation during ablation of porcine skin with erbium:YAG laser vs a novel picosecond infrared laser. *JAMA Otolaryngol Head Neck Surg*. 2013;139(8):828-833.

- [21] Jowett N, Wollmer W, Reimer R, et al. Bone ablation without thermal or acoustic mechanical injury via a novel picosecond infrared laser (PIRL). *Otolaryngol Head Neck Surg*. 2014;150(3):385-393.
- [22] Geronemus RG, Kazlouskaya V, Bae YS, et al. Diffractive lens array with picosecond laser for facial acne scarring: follow up and histology. Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [23] Weiss RA, Weiss MA, Beasley K. Picosecond laser for reduction of wrinkles. Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [24] McDaniel D. Treatment of facial photodamage and rhytides using a picosecond pulsed alexandrite laser and specifically designed focus Optico Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [25] Kilmer SL, Custis, T. Single vs repeat exposure tattoo removal during single sessions with pico-second pulse duration laser technology. Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [26] Tanghetti E, Tanghetti M. Dose optimization with a picosecond 755nm alexandrite laser for tattoo removal. Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [27] Alabulrazzaq H, Bae YSC, Brauer JA, Geronemus RG. Rapid removal of red tattoo pigment with a novel 532nm Nd:YAG picosecond laser. Paper presented at: 34th Annual Scientific Conference of the American Society for Laser Medicine & Surgery; April 2014; Phoenix, AZ.
- [28] Diehl J, Chesnut C, Lask G. Treatment of nevus of ota with picosecond alexandrite laser and review of the dermal melanocytoses and treatment options. Treatment of Nevus of Ota with picosecond Alexandrite Laser and review of the dermal melanocytoses and treatment options. American Society for Laser Medicine & Surgery 2014;Annual Conference 2014, Phoenix, AZ.
- [29] Tanghetti E, Tanghetti M. A clinical and histological study of skin treated with a pico-second alexandrite laser comparing a uniform treatment spot and a spatially modulated spot. Treatment of Nevus of Ota with picosecond Alexandrite Laser and review of the dermal melanocytoses and treatment options. American Society for Laser Medicine & Surgery 2014;Annual Conference 2014, Phoenix, AZ.

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- [1] Ho C, Nguyen Q, Lowe NJ 等人著, Laser Resurfacing in Pigmented Skin. *Dermatology Surg* 1995;21:1035-1037.
- [2] Nelson S 著, Lasers and Laser-Tissue Interactions in Dermatology. Principles and Practices in Cutaneous Laser Surgery. Kauvar A, Hruza G. 2005:66.
- [3] Anderson RR, Margolis RJ, Dover J 等人著, Selective Photothermolysis of cutaneous pigmentation by Q-Switched Nd:YAG laser pulsed at 1064 nm, 532 nm, and 355 nm. *J Invest Dermatol* 1989;93:28-32.
- [4] Dover J, Arndt K, Ort R 著, Lasers in the treatment of pigmented lesions. Principles and Practices in Cutaneous Laser Surgery. Kauvar A, Hruza G. 2005:66.
- [5] Saedi N, Metelitsa A, Petrell K 等人著, Treatment of Tattoos With a Picosecond Alexandrite Laser: A Prospective Trial. *Arch Dermatol*. 2012;148(12):1360-1363.
- [6] Olson RL, Gaylor J, Everett MA 著, Skin color, melanin, and erythema. *Arch Dermatology* 1988; 124:869.
- [7] Brauer JA, Kazlouskaya V, Alabdulrazzaq H 等人著, Use of a Picosecond Pulse Duration Laser With Specialized Optic for Treatment of Facial Acne Scarring. *JAMA Dermatol*. Published online November 19, 2014.

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- [1] Rabinovitz HS, Barnhill RL. Benign melanocytic neoplasms. In: Bologna JL, Jorizzo JL, Schaffer JV, editors. Dermatology. 3rd ed. London: Elsevier Saunders; 2012; pp. 1858-9.
- [2] Anderson RR, Parrish JA. Selective photothermolysis: precise microsurgery by selective absorption of pulsed radiation. Science 1983; 220:524-7.
- [3] Alster TS, Williams CM. Treatment of nevus of Ota by the Q-switched alexandrite laser. Dermatol Surg 1995;21:592-6.
- [4] Felton SJ, Al-Niami F, Ferguson JE, Madan V. Our perspective of the treatment of nevus of Ota with 1,064-, 755- and 532-nm wavelength lasers. Lasers Med Sci 2014;29:1745-9.
- [5] Ross V, Naseef G, Lin G, et al. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. Arch Dermatol 1998;134:167-71.

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- [1] Izikson L, Farinelli W, Sakamoto F, Tannous Z, Anderson RR. Safety and effectiveness of black tattoo clearance in a pig model after a single treatment with a novel 758nm 500 picosecond laser: A pilot study. *Lasers Surg Med* 2010; 42(7): 640-646.
- [2] Ross V, Naseef G, Lin G, Kelly M, Michaud N, Flotte TJ, Raythen J, Anderson RR. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. *Arch Dermatol* 1998;134:167-171.
- [3] Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: A prospective trial. *Arch Dermatol* 2012;148:1360-1363.
- [4] Brauer JA, Reddy KK, Anolik R, Weiss ET, Karen JK, Hale EK, Brightman LA, Bernstein L, Geronemus RG. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. *Arch Dermatol* 2012;148:820-823.
- [5] Tanghetti EA, Tanghetti M. Dose optimization with a picosecond 755nm alexandrite laser for tattoo removal. 34th American Society for Laser Medicine and Surgery (ASLMS) annual conference. Phoenix Arizona. 2014.
- [6] Kilmer S, Custis T. Single vs repeat exposure tattoo removal during single sessions with picosecond pulse duration laser technology. 34th American Society for Laser Medicine and Surgery (ASLMS) annual conference Phoenix Arizona. 2014.
- [7] Chesnut C, Diehl J, Lask G. Treatment of Nevus of Ota with a picosecond 755-nm alexandrite laser. *Dermatol Surg* 2015; 41(4):508-510.
- [8] Kono T, Nozaki M, Chan HH, Mikashima Y. A retrospective study looking at the long-term complications of Q-switched ruby laser in the treatment of nevus of Ota. *Lasers Surg Med* 2001;29:156-159.
- [9] Chan HH, Ying SY, Ho WS, Kono T, King WW. An in vivo trial comparing the clinical efficacy and complications of Q-switched 755 nm alexandrite and Q-switched 1064 nm Nd:YAG lasers in the treatment of nevus of Ota. *Dermatol Surg* 2000;26:919-922.
- [10] Changzheng H, Sami LS, Yan L. Factors affecting response and laser sessions and complications in Nevus of Ota treated by Q-switched alexandrite laser: A retrospective study. *G Ital Dermatol Venereol* 2015.
- [11] Felton SJ, Al-Niami F, Ferguson JE, Madan V. Our perspective of the treatment of nevus of Ota with 1,064-, 755- and 532-nm wavelength lasers. *Lasers Med Sci* 2014;29:1745-1749.
- [12] Kilmer SL, Wheeland RG, Goldberg DJ, Anderson RR. Treatment of epidermal pigmented lesions with the frequency-doubled Q-switched Nd:YAG laser. A controlled, single-impact, dose-response, multicenter trial. *Arch Dermatol* 1994;130:1515.
- [13] Kagami S, Asahina A, Watanabe R, Mimura Y, Shirai A, Hattori N, Watanabe T, Tamaki K. Treatment of 153 Japanese patients with Q-switched alexandrite laser. *Lasers Med Sci* 2007;22:159-163.
- [14] Wang Y, Qian H, Lu Z. Treatment of café-au-lait macules in Chinese patients with a Q-switched 755-nm alexandrite laser. *J Dermatol Treat* 2012;23(6):431-436.
- [15] Kopf AW, Weidman AI. Nevus of Ota. *Arch Dermatol* 1962;85:195-208.
- [16] Hori Y, Oohara K, Niimura M, Kukita A. Electron microscopy: Ultrastructural observations of the extracellular sheath of dermal melanocytes in the nevus of Ota. *Am J Dermatopathol* 1982;4:245-251.
- [17] Wang CC, Sue YM, Yang CH, Chen CK. A comparison of Q-switched alexandrite laser and intense pulsed light for the treatment of freckles and lentigines in Asian persons: a randomised, physician-blinded, split-face comparative trial. *J Am Acad Dermatol* 2006;54:804-810.
- [18] Murphy MJ, Huang MY. Q-switched ruby laser treatment of benign pigmented lesions in Chinese skin. *Ann Acad Med Singapore* 1994;23:60-66.
- [19] Chan HH, Fung WK, Ying SY, Kono T. An in vivo trial comparing the use of different types of 532 nm Nd:YAG lasers in the treatment of facial lentigines in oriental patients. *Dermatol Surg* 2000;26(8):743-749.
- [20] Choi JE, Kim JW, Seo SH, Son SW, Ahn HH, Kye YC. Treatment of Becker's nevi with a long-pulsed alexandrite laser. *Dermatol Surg* 2009;35:1105-1108.
- [21] Trelles MA, Allones I, Moreno-Arias GA, Velez M. Becker's naevus: A comparative study between erbium:YAG and Q-switched neodymium:YAG; clinical and histopathological findings. *Br J Dermatol* 2005;152:308-313.
- [22] Kunachak S, Leelaudomlapi P, Sirikulchayanonta V. Q-switched ruby laser therapy of acquired bilateral nevus of Ota-like macules. *Dermatol Surg* 1999;25:938-934.
- [23] Lam AY, Wong DS, Lam LK, Ho WS, Chan HH. A retrospective study on the efficacy and complications of Q-switched alexandrite laser in the treatment of acquired bilateral nevus of Ota-like macules. *Dermatol Surg* 2001;27: 937-941.
- [24] Polnikom N, Tanrattanakom S, Goldberg DJ. Treatment of Hori's nevus with the Q-switched Nd:YAG laser. *Dermatol Surg* 2000;26:477-480.
- [25] Wang CC, Sue YM, Yang CH, Chen CK. A comparison of Q-switched alexandrite laser and intense pulsed light for the treatment of freckles and lentigines in Asian persons: A randomised, physician-blinded, split-face comparative trial. *J Am Acad Dermatol* 2006;54:804-810.
- [26] Murphy MJ, Huang MY. Q-switched ruby laser treatment of benign pigmented lesions in Chinese skin. *Ann Acad Med Singapore* 1994;23:60-66.
- [27] Chan HH, Fung WKK, Ying SY, Kono T. An in vivo trial comparing the use of different types of 532-nm Nd:YAG lasers in the treatment of facial lentigines in oriental patients. *Dermatol Surg* 2000;26:743-749.
- [28] Negishi K, Tezuka Y, Kushikata N, Wakamatsu S. Photorejuvenation for Asian skin by intense pulsed light. *Dermatol Surg* 2001;27:627-631.
- [29] Ho SG, Chan PY, Yeung CK, Shek SY, Kono T, Chan HH. A retrospective analysis of the management of freckles and lentigines using four different pigment lasers on Asian skin. *J Cosmet Laser Ther* 2012;14:74-80.

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- Alexiades-Armenakas M. Rhytides, laxity, and photoaging treated with a combination of radio frequency, diode laser, and pulsed light and assessed with a comprehensive grading scale. *J Drugs Dermatol* 2006;5(8):731-738.
- Alexiades-Armenakas MR, Dover JS, Arndt KA. Fractional laser skin resurfacing. *J Drugs Dermatol* 2012;11(11):1274-1287.
- Rabe JH, Mamelak AJ, McElgunn PJ, Morison WL, Sauder DN. Photoaging: Mechanisms and repair. *J Am Acad Dermatol* 2006;55(1):1-19.
- Gold MH, Sensing W, Biron J. Fractional Q-switched 1,064-nm laser for the treatment of photoaged-photodamaged skin. *J Cosmet Laser Ther* 2014;16(2):69-76.
- Huo MH, Wang YQ, Yang X. Split-face comparison of intense pulsed light and nonablative 1,064-nm Q-switched laser in skin rejuvenation. *Dermatol Surg* 2011;37(1):52-57.
- Luebberding S, Alexiades-Armenakas MR. Fractional, nonablative Q-switched 1,064-nm neodymium YAG laser to rejuvenate photoaged skin: A pilot case series. *J Drugs Dermatol* 2012;11(11):1300-1304.
- Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. treatment of tattoos with a picosecond alexandrite laser: A prospective trial. *Arch Dermatol* 2012;148(12):1360-1363.
- Chesnut C, Diehl J, Lask G. Treatment of Nevus of Ota with a Picosecond 755-nm Alexandrite laser. *Dermatol Surg* 2015;41(4):508-510.
- Brauer JA, Kazlouskaya V, Alabdulrazzaq H, Bae YS, Bernstein LJ, Anolik R, Heller PA, Geronemus RG. Use of a picosecond pulse duration laser with specialized optic for treatment of facial acne scarring. *JAMA Dermatol* 2015;151(3):278-284.
- Alexiades-Armenakas M. A quantitative and comprehensive grading scale for rhytides, laxity, and photoaging. *J Drugs Dermatol* 2006;5(8):808-809.
- Wu DC, Friedmann DP, Fabi SG, Goldman MP, Fitzpatrick RE. Comparison of intense pulsed light with 1,927-nm fractionated thulium fiber laser for the rejuvenation of the chest. *Dermatol Surg* 2014;40(2):129-133.

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- [1] Halder RM, Holmes YC, Bridgemen-Shah S, Kligman AM. A clinicohistopathologic study of acne vulgaris in black females. *J Invest Dermatol* 1996;106:495A.
- [2] Alexis AF, Sergay AB, Taylor SC. Common dermatologic disorders in skin of color: A comparative practice survey. *Cutis* 2007;80:387-394.
- [3] United States Census Bureau. Census Explorer: Population Estimates Edition: Total Population. <http://www.census.gov/censusexplorer/censusexplorer-popest.html>. Washington, DC: Populations projections program, population division, U.S. Census Bureau.
- [4] Projections of the resident population by race, hispanic origin, and nativity: Middle series, 2050 to 2070. <https://www.census.gov/population/projections/files/natproj/summary/np-t5-g.pdf>. Washington, DC: Populations projections program, population division, U.S. Census Bureau.
- [5] Ortonne JP, Pandya AG, Lui H, Hexsel D. Treatment of solar lentigines. *J Am Acad Dermatol* 2006;54(5 Suppl 2):S262-271.
- [6] Kunachak S, Kunacharkr S, Sirikulchayanonta V, Leelaudomniti P. Dermabrasion is an effective treatment for acquired bilateral nevus of Ota-like macules. *Dermatol Surg* 1996;22:559-562.
- [7] Anderson RR, Parish JA. Selective photothermolysis: Precise microsurgery by selective absorption of pulsed radiation. *Science* 1983;220(4596):524-527.
- [8] Geronemus RG. Fractional photothermolysis: Current and future applications. *Lasers Surg Med* 2006;38(3):169-176.
- [9] Anderson RR, Margolis RJ, Watanabe S, Flotte T, Hruza GJ, Dover JS. Selective photothermolysis of cutaneous pigmentation by Q-switched Nd:YAG laser pulses at 1064, 532, and 355nm. *J Invest Dermatol* 1989;93(1):28-32.
- [10] Saedi N, Metelitsa A, Petreel K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: A prospective trial. *Arch Dermatol* 2012;148:1360-1363.
- [11] Brauer JA, Reddy KK, Anolik R, Weiss ET, Karen JK, Hale EK, Brightman LA, Bernstein L, Geronemus RG. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. *Arch Dermatol* 2012;148:820-823.
- [12] Grevelink JM, van Leeuwen RL, Anderson RR, Byers HR. Clinical and histological responses of congenital melanocytic nevi after single treatment with Q-switched lasers. *Arch Dermatol* 1997;133:349-353.
- [13] Sadigha A, Saatee S, Muhagheh-Zahed G. Efficacy and adverse effects of Q-switched ruby laser on solar lentigines: A prospective study of 91 patients with Fitzpatrick skin type II, III, and IV. *Arch Dermatol* 1997;133:349-353.
- [14] Chestnut C, Diehl J, Lask G. Treatment of nevus of Ota with a picosecond 755-nm alexandrite laser. *Dermatol Surg* 2015;41: 508-510.
- [15] Polla L, Margolis R, Dover JS, Whitaker D, Murphy GF, Jacques SL, Anderson SS. Melanosomes are a primary target of Q-switched ruby laser irradiation in guinea pig skin. *J Invest Dermatol* 1987;89(3):281-286.
- [16] Ross V, Naseef G, Lin C, Kelly M, Michaud N, Flotte T, Raythen J, Anderson R. Comparison of responses of tattoos to picosecond and nanosecond Q-switched Neodymium:YAG lasers. *Arch Dermatol* 1998;134:167-171.
- [17] DD Ho, R Lond on, GB Zimmerman, DA. Young Laser-tattoo removal: A study of the mechanism and the optimal treatment strategy via computer stimulations. *Lasers Surg Med* 2002;30(5):389-397.

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- 2013;15(2):70-73.
- [11] Kim BJ, Lee DH, Kim MN, et al. Fractional photothermolysis for the treatment of striae distensae in Asian skin. *Am J Clin Dermatol.* 2008;9(1):33-37.
- [12] Katz TM, Goldberg LH, Friedman PM. Nonablative fractional photothermolysis for the treatment of striae rubra. *Dermatol Surg.* 2009;35(9):1430-1433.
- [13] Yang YJ, Lee G-Y. Treatment of striae distensae with nonablative fractional laser versus ablative CO<sub>2</sub> fractional laser: a randomized controlled trial. *Ann Dermatol.* 2011;23(4):481-489.
- [14] Lee HM, Haw S, Kim JK, Chang SE, Lee MW. Split-face study using a 1,927-nm thulium fiber fractional laser to treat photoaging and melasma in Asian skin. *Dermatol Surg.* 2013;39(6):879-888.
- [15] Toulaki A, Galimberti MG, Pellacani G, Bencini PL. Combination of fractional erbium-glass laser and topical therapy in melasma resistant to triple-combination cream. *J Dermatolog Treat.* 2014;25(3):218-222.
- [16] Chan NP, Ho SG, Yeung CK, Shek SY, Chan HH. Fractional ablative carbon dioxide laser resurfacing for skin rejuvenation and acne scars in Asians. *Lasers Surg Med.* 2010;42(9): 615-623.
- [17] Alexis AF. Fractional laser resurfacing of acne scarring in patients with Fitzpatrick skin types IV-VI. *J Drugs Dermatol.* 2011;10(12 Suppl):s6-s7.
- [18] Clark CM, Silverberg JJ, Alexis AF. A retrospective chart review to assess the safety of nonablative fractional laser resurfacing in Fitzpatrick skin types IV to VI. *J Drugs Dermatol.* 2013;12(4): 428-431.
- [19] Wang C-C, Sue Y-M, Yang C-H, Chen C-K. A comparison of Q-switched alexandrite laser and intense pulsed light for the treatment of freckles and lentigines in Asian persons: a randomized, physician-blinded, split-face comparative trial. *J Am Acad Dermatol.* 2006;54(5):804-810.
- [20] Watanabe S, Takahashi H. Treatment of nevus of Ota with the Q-switched ruby laser. *N Engl J Med.* 1994;331(26):1745-1750.
- [21] Polnikorn N, Tanrattanakorn S, Goldberg DJ. Treatment of Hori's nevus with the Q-switched Nd:YAG laser. *Dermatol Surg.* 2000;26(5):477-480.
- [22] Ho SG, Chan HH. The Asian dermatologic patient: review of common pigmentary disorders and cutaneous diseases. *Am J Clin Dermatol.* 2009;10(3):153-168.
- [23] Sadighha A, Saatei S, Muhagheh-Zahed G. Efficacy and adverse effects of Q-switched ruby laser on solar lentigines: a prospective study of 91 patients with Fitzpatrick skin type II, III, and IV. *Dermatol Surg.* 2008;34(11):1465-1468.
- [24] Chan HH, Fung WK, Ying SY, Kono T. An in vivo trial comparing the use of different types of 532 nm Nd:YAG lasers in the treatment of facial lentigines in Oriental patients. *Dermatol Surg.* 2000;26(8):743-749.
- [25] Murphy MJ, Huang MY. Q-switched ruby laser treatment of benign pigmented lesions in Chinese skin. *Ann Acad Med Singap.* 1994;23(1):60-66.
- [26] Ross V, Naseef G, Lin G, et al. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. *Arch Dermatol.* 1998;134(2):167-171.
- [27] Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: a prospective trial. *Arch Dermatol.* 2012;148(12):1360-1363.
- [28] Brauer JA, Reddy KK, Anolik R, et al. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. *Arch Dermatol.* 2012;148(7):820-823.
- [29] Fabi SG, Metelitsa AI. Future directions in cutaneous laser surgery. *Dermatol Clin.* 2014;32(1):61-69.
- [30] The American Society for Aesthetic Plastic Surgery. ASAPS 1999 Statistics on Cosmetic Surgery. Available from: URL: <http://www.surgery.org/sites/default/files/ASAPS1999Stats.pdf>. Accessed September 21, 2015.
- [31] Cosmetic Surgery National Data Bank. Statistics 2013. *Aesthet Surg J.* 2014;34(Suppl):1S-22S.
- [32] Brauer JA, Kazlouskaya V, Alabdulrazzaq H, et al. Use of a picosecond pulse duration laser with specialized optic for treatment of facial acne scarring. *JAMA Dermatol.* 2015;151(3): 278-284.
- [33] Tanghetti E, Tanghetti M. A clinical and histologic study of skin treated with a picosecond alexandrite laser comparing a uniform treatment. *Lasers Surg Med.* 2014;46(S25):28.
- [34] Tanghetti E. Characterization of the histologic changes in the skin from the treatment with the 755nm picosecond alexandrite laser. *Lasers Surg Med.* 2015;47(S26):24.
- [35] Sommer S, Seukeran DC, Sheehan-Dare RA. Efficacy of pulsed dye laser treatment of port wine stain malformations of the lower limb. *Br J Dermatol.* 2003;149(4):770-775.
- [36] Grimes PE. Management of hyperpigmentation in darker racial ethnic groups. *Semin Cutan Med Surg.* 2009;28(2):77-85.
- [37] Cheyasak N, Manuskiatti W, Maneeprasopchoke P, Wanitphakdeedecha R. Topical corticosteroids minimize the risk of postinflammatory hyperpigmentation after ablative fractional CO<sub>2</sub> laser resurfacing in Asians. *Acta Derm Venereol.* 2015;95(2):201-205.

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- [1] ASDS [webpage on the internet]. Laser Therapy for Unwanted Tattoos. Chicago: ASDS Public Resources; c1971-2016 [updated 2013; cited August 16, 2015]. Available from: [https://www.asds.net/\\_PublicResources.aspx?id=6073](https://www.asds.net/_PublicResources.aspx?id=6073). Accessed August 16, 2015.
- [2] Kent KM, Graber EM. Laser tattoo removal: a review. *Dermatol Surg.* 2012;38(1):1-13.
- [3] Choudhary S, Elsaie ML, Leiva A, Nouri K. Lasers for tattoo removal: a review. *Lasers Med Sci.* 2010;25(5):619-627.
- [4] Ho SG, Goh CL. Laser tattoo removal: a clinical update. *J Cutan Aesthet Surg.* 2015;8(1):9-15.
- [5] Anderson RR, Parrish JA. Selective photothermolysis: precise micro-surgery by selective absorption of pulsed radiation. *Science.* 1983;220(4596):524-527.
- [6] Ho DD, London R, Zimmerman GB, Young DA. Laser-tattoo removal—a study of the mechanism and the optimal treatment strategy via computer simulations. *Lasers Surg Med.* 2002;30(5):389-397.
- [7] Ross V, Naseef G, Lin G, et al. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. *Arch Dermatol.* 1998;134(2):167-171.
- [8] Herd RM, Alora MB, Smoller B, Arndt KA, Dover JS. A clinical and histologic prospective controlled comparative study of the picosecond titanium:sapphire (795 nm) laser versus the Q-switched alexandrite (752 nm) laser for removing tattoo pigment. *J Am Acad Dermatol.* 1999;40(4):603-606.
- [9] PicoSure.com [homepage on the Internet]. Westford: PicoSure ; Specifications, Cynosure Inc.; c2014 [updated 2014; cited August 16, 2015]. Available from: <http://www.cynosure.com/product/picosure/>. Accessed August 16, 2015.
- [10] Izikson L, Farinelli W, Sakamoto F, Tannous Z, Anderson RR. Safety and effectiveness of black tattoo clearance in a pig model after a single treatment with a novel 758 nm 500 picosecond laser: a pilot study. *Lasers Surg Med.* 2010;42(7):640-646.
- [11] Brauer JA, Reddy KK, Anolik R, et al. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. *Arch Dermatol.* 2012;148(7):820-823.
- [12] Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: a prospective trial. *Arch Dermatol.* 2012;148(12):1360-1363.
- [13] Au S, Liolios AM, Goldman MP. Analysis of incidence of bulla formation after tattoo treatment using the combination of the picosecond Alexandrite laser and fractionated CO<sub>2</sub> ablation. *Dermatol Surg.* 2015;41(2):242-245.
- [14] Alabdulrazzaq H, Brauer JA, Bae YS, Geronemus RG. Clearance of yellow tattoo ink with a novel 532-nm picosecond laser. *Lasers Surg Med.* 2015;47(4):285-288.
- [15] Bernstein EF, Schomacker KT, Basilavecchio LD, Plugis JM, Bhawalkar JD. A novel dual-wavelength, Nd:YAG, picosecond-domain laser safely and effectively removes multicolor tattoos. *Lasers Surg Med.* Epub 2015 Jul 14.
- [16] Anderson RR, Donelan MB, Hivnor C, et al. Laser treatment of traumatic scars with an emphasis on ablative fractional laser resurfacing: consensus report. *JAMA Dermatol.* 2014;150(2):187-193.

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- [1] Goldman L, Rockwell RJ, Meyer R, Otten R, Wilson RG, Kitzmiller KW. Laser treatment of tattoos: a preliminary survey of three year 's clinical experience. JAMA. 1967;201(11):841-844.
- [2] Taylor CR, Gange RW, Dover JS, et al. Treatment of tattoos by Q-switched ruby laser: a dose-response study. Arch Dermatol. 1990;126(7):893-899.
- [3] Anderson RR, Parrish JA. Selective photothermolysis: precise microsurgery by selective absorption of pulsed radiation. Science. 1983;220(4596):524-527.
- [4] Ross V, Naseef G, Lin G, et al. Comparison of responses of tattoos to picosecond and nanosecond Q-switched neodymium: YAG lasers. Arch Dermatol. 1998;134 (2):167-171.
- [5] Herd RM, Alora MB, Smoller B, Arndt KA, Dover JS. A clinical and histologic prospective controlled comparative study of the picosecond titanium:sapphire (795 nm) laser versus the Q-switched alexandrite (752 nm) laser for removing tattoo pigment. J Am Acad Dermatol. 1999;40(4):603-606.
- [6] Izikson L, Farinelli W, Sakamoto F, Tannous Z, Anderson RR. Safety and effectiveness of black tattoo clearance in a pig model after a single treatment with a novel 758 nm 500 picosecond laser: a pilot study. Lasers Surg Med. 2010;42 (7):640-646.
- [7] Fitzpatrick RE, Goldman MP. Tattoo removal using the alexandrite laser. Arch Dermatol. 1994;130(12):1508-1514.
- [8] Alster TS. Q-switched alexandrite laser treatment (755 nm) of professional and amateur tattoos. J Am Acad Dermatol. 1995;33(1):69-73.

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- [1] Goldman L, Rockwell RJ, Meyer R, Otten R, Wilson RG, Kitzmiller KW. Laser treatment of tattoos: a preliminary survey of three years' clinical experience. *JAMA*. 1967;201(11):841-844.
- [2] Reid R, Muller S. Tattoo removal by CO laser dermabrasion. *Plast Reconstr Surg*. 1980;65(6):717-728.
- [3] Apfelberg DB, Maser MR, Lash H. Argon laser treatment of decorative tattoos. *Br J Plast Surg*. 1979;32(2):141-144.
- [4] Anderson RR, Parrish JA. Microvasculature can be selectively damaged using dye lasers: a basic theory and experimental evidence in human skin. *Lasers Surg Med*. 1981;1(3):263-276.
- [5] Anderson RR, Parrish JA. Selective photothermolysis: precise microsurgery by selective absorption of pulsed radiation. *Science*. 1983;220(4596):524-527.
- [6] Reid WH, McLeod PJ, Ritchie A, Ferguson-Pell M. Q-switched ruby laser treatment of black tattoos. *Br J Plast Surg*. 1983;36(4):455-459.
- [7] Graber E, Iyengar V, Rohrer T, et al. Laser treatment of tattoos and pigmented lesions. In: Robinson JK, Hanke CW, Siegel DM, Fratila A, eds. *Surgery of the Skin: Procedural Dermatology*. 2nd ed. China: Mosby; 2010:537-548.
- [8] Ara G, Anderson RR, Mandel KG, Ottesen M, Oseroff AR. Irradiation of pigmented melanoma cells with high intensity pulsed radiation generates acoustic waves and kills cells. *Lasers Surg Med*. 1990;10(1):52-59.
- [9] Kent KM, Graber EM. Laser tattoo removal: a review. *Dermatol Surg*. 2012;38(1): 1-13.
- [10] Taylor CR, Anderson RR, Gange RW, Michaud NA, Flotte TJ. Light and electron microscopic analysis of tattoos treated by Q-switched ruby laser. *J Invest Dermatol*. 1991;97(1):131-136.
- [11] Leuenberger ML, Mulas MW, Hata TR, Goldman MP, Fitzpatrick RE, Grevelink JM. Comparison of the Q-switched alexandrite, Nd:YAG, and ruby lasers in treating blue-black tattoos. *Dermatol Surg*. 1999;25(1):10-14.
- [12] Lin T, Jia G, Rong H, Li J, Zhou Z. Comparison of a single treatment with Q-switched ruby laser and Q-switched Nd:YAG laser in removing black-blue Chinese tattoos. *J Cosmet Laser Ther*. 2009;11(4):236-239.
- [13] Zelickson BD, Mehregan DA, Zarrin AA, et al. Clinical, histologic, and ultrastructural evaluation of tattoos treated with three laser systems. *Lasers Surg Med*. 1994;15(4):364-372.

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转载: CDR E. Victor Ross, USN, Department of Clinical Research, Naval Medical Center, 34800 Bob Wilson Dr, San Diego, CA 92134-5000.

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- [1] Vitkin I, Wilson B, Anderson R, Prahl S. Pulsed photothermal radiometry in optically transparent media containing discrete optical absorbers. *PhysMedBiol*. 1994;39:1-24.
- [2] Kilmer SL, Lee MS, Grevelink JM, Flotte TJ, Anderson RR. The Q-switched Nd: YAG laser effectively treats tattoos. *Arch Dermatol*. 1993;129:971-978.
- [3] Kilmer SL, Anderson RR. Clinical use of the Q-switched ruby and the Q-switched Nd:YAG (1064 nm and 532 nm) lasers for treatment of tattoos. *JDermatolSurgOncol*. 1993;19:330-338.
- [4] Chen H, Diebold G. Chemical generation of acoustic waves: a giant photoacoustic effect. *Science*. 1995;270:963-966.
- [5] Anderson R. Laser-tissue interactions. In: Goldman M, Fitzpatrick R, eds. *CutaneousLaserSurgery*. St Louis, Mo: Mosby Year Book Inc; 1994:1-9.
- [6] Taylor CR, Anderson RR, Gange RW, Michaud NA, Flotte TJ. Light and electron microscopic analysis of tattoos treated by Q-switched ruby laser. *JInvestDermatol*. 1991;97:131-136.
- [7] Taylor CR, Gange RW, Dover JS, et al. Treatment of tattoos by Q-switched ruby laser: a dose-response study. *ArchDermatol*. 1990;126:893-899.
- [8] Zelickson BD, Mehregan DA, Zarrin AA, et al. Clinical, histologic, and ultrastructural evaluation of tattoos treated with three laser systems. *LasersSurgMed*. 1994;15:364-372.
- [9] Watanabe S, Anderson R, Brorson S, Daliakas G, Fujimoto J, Flotte T. Comparative studies of femtosecond to microsecond laser pulses on selective pigmented cell injury in skin. *PhotochemPhotobiol*. 1991;53:757-762.
- [10] Carslaw H, Jaeger J. *HeatConductioninSolids*. 2nd ed. Oxford, England: Clarendon Press; 1986.
- [11] van Leewen T, Jansen E, Motamedi M, Borst C, Welch A. Pulse laser ablation of soft tissue. In: Welch A, van Gemert M, eds. *Optical-ThermalResponseofLaserIrradiatedTissue*. New York, NY: Plenum Publishing Corp; 1995:712-717.
- [12] Doukas A, Flotte T. Physical characteristics and biological effects of laser-induced stress waves. *UltrasoundMedBiol*. 1996;22:1-9.

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- [1] Brauer JA, Reddy KK, Anolik R, Weiss ET, Karen JK, Hale EK, Brightman LA, Bernstein L, Geronemus RG. Successful and rapid treatment of blue and green tattoo pigment with a novel picosecond laser. Arch Dermatol 2012;148(7):820-823.
- [2] Saedi N, Metelitsa A, Petrell K, Arndt KA, Dover JS. Treatment of tattoos with a picosecond alexandrite laser: A prospective trial. Arch Dermatol 2012;148(12):1360-1363.
- [3] Kossida T, Rigopoulos D, Katsambas A, Anderson RR. Optimal tattoo removal in a single laser session based on the method of repeated exposures. J Am Acad Dermatol 2012;66(2):271-277.
- [4] Reddy KK, Brauer JA, Anolik R, Bernstein L, Brightman L, Hale E, Karen J, Weiss E, Geronemus RG. Topical per urodecalin resolves immediate whitening reactions and allows rapideffectivemultiplepasstreatmentoftattoos. Lasers Surg Med 2013;45(2):76-80.
- [5] Bencini PL, Cazzaniga S, Tourlaki A, Galimberti MG, Naldi L. Removal of tattoos by q-switched laser: variables influencing outcome and sequelae in a large cohort of treated patients. Arch Dermatol 2012;148(12):1364-1369.
- [6] Kirby W, Chen CL, Desai A, Desai T. Causes and recommendations for unanticipated ink retention following tattoo removal treatment. J Clin Aesthet Dermatol 2013;6(7):27-31.
- [7] Ferguson JE, August PJ. Evaluation of the Nd:YAG laser for the treatment of amateur and professional tattoos. Br J Dermatol 1996;135(4):586-591.
- [8] Ho DD, London R, Zimmerman GB, Young DA. Laser-tattoo removal-a study of the mechanism and the optimal treatment strategy via computer simulations. Lasers Surg Med 2002;30(5):389-397.
- [9] Gomez C, Martin V, Sastre R, Costela A, Garca-Moreno I. In vitro and in vivo laser treatments of tattoos: High efficiency and low sequelae. Arch Dermatol 2010;146(1):39-45.
- [10] Beute TC, Miller CH, Timko AL, Ross EV. In vitro spectral analysis of tattoo pigments. Dermatol Surg 2008;34(4):508-515. discussion 515-6.

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